

Product Name: Electric Sieving Apparatus with JJSJ Grain Sieve

Characteristic

Electrically screening instrument uses the table model structure, the body sieves by three planes struts, makes the plane rotation through the eccentric link motion gear the automatic device to sieve, has the structure compact reasonable, the revolution steady, the performance is stable merit and so on ease of operation, widely adapts the grain, the fuel oils screening test.

Sieve is suitable to measure the impurity percentage in grain sample or the bran percentage in rice. It is an inspecting instrument to judge the samples quality. Sieve is suitable to measure the insect pest's percentage in grain and oil plant, and determine insect density and judge grains grade. The sieves are made on aluminum. It can be used with our electric screen meter.



Technical Specification (Electric Sieving Apparatus)

Max sieving quantity range: 500g

Sieving frame level: 3 layers

Sieving range: 100mm

Rotation speed: 115±5r/min

Sieve diameter: 220mm

Power: 60W AC220V 50Hz

Weight: 13kg

Package Size: 42*45*31CM

Technical Specification (Grain Sieve):

Sieve Layer: 12

Diameter: 220mm

Height: 50mm

Sieve opening dimensions(mm) : Φ 1.0 ,1.5, 2.0, 2.5, 3.0, 3.5, 4.0, 4.5, 5.0, 6.0, 1.7× 20, 12

Packing: standard exporting package

Related products

Automatic Seed Counter/Grain Moisture Meter(with volume-weight)/Grain Moisture Meter/Seed Blower